Thermo Scientific

mySPIN 6 Microcentrifuge

Instruction Manual

50145670-e • 01/ 2025

Instruction Manual	2
Manual de Instrucciones	12
Manual d'Instructions.....	22
Bedienungsanleitung.....	32
Manuale de Istruzioni	42

Safety Precautions

The Thermo Scientific™ mySPIN 6 Microcentrifuge is a personal mini-centrifuge designed for laboratories in the bioscience, medical, and chemistry fields. It may be used with biohazardous and/or hazardous materials. This instruction manual cannot address all safety hazards. It is the responsibility of the user to consult and observe all health and safety precautions and to assess the instrument's suitability to the task.

CAUTION

Before starting the centrifuge for the first time, please read this instruction manual carefully. This manual contains important information, safeguards, and operating instructions.

1. Never use this product in any manner inconsistent with these instructions.
2. This product is intended to be used indoors only.
3. This product is intended for separating aqueous solutions in approved test tubes compatible to the RCF of this unit.
4. This product is not intended to centrifuge very dense materials. Do not use materials with a density greater than 1.2g/cm³ (1.2g/mL).
5. Use only the power adapter delivered with the product. Use of other power adapters will void warranty.
6. Do not put hands in the rotor area unless the rotor is completely stopped.
7. Only use the originally supplied rotors, AC adapters, tube adapters and any other spare parts delivered with the unit. These parts have been designed to safely meet the maximum RCF that the centrifuge is able to deliver. **DO NOT USE ROTORS OR COMPONENTS FROM ANY OTHER MANUFACTURER.**
8. Always make sure the rotors are securely seated.
9. Never use this product in any manner inconsistent with these instructions.
10. Observe general laboratory safety precautions and regulations.
11. Do not attempt to operate with the cover removed or open.
12. The rotor and lid must always be securely fastened when in operation.
13. Do not use damaged rotors.
14. Do not move unit while the rotor is spinning.
15. Do not attempt to stop the rotor while the unit is running. Doing so may cause injury or the unit to fail and will void the warranty.
16. Do not immerse unit in water. Injury due to shock or fire may result.
17. The rotor must be loaded symmetrically. Operating this unit with an unbalanced rotor will cause severe vibration and damage, which is not covered by warranty.
18. Do not fill tubes while inserted in rotor. Spilled liquid may harm unit.

19. Close tube lids before starting centrifuge. Open lids can be ripped off during centrifugation and damage the unit. Open lids can cause fluid to be dispersed inside the centrifuge.
20. Do not insert non-standard tubes or other unfamiliar articles into the rotor.
21. Follow all precautions from the Material Safety Data Sheet (MSDS) for any reagent you use with this instrument.
22. Do not use solvent or flammable liquids near this or other electrical equipment.
23. Do not centrifuge materials that in combination or alone are flammable, explosive, or corrosive, or may cause any other hazard.
24. Do not operate in a hazardous or flammable environment.
25. Do not use any solvent on the unit that may attack plastic or rubber or cause cracks in the rotor.
26. Ensure the rotor is protected from corrosion and mechanical damage. The rotor must be cleaned with a pH-neutral cleaning liquid.
27. Before using any cleaning or decontamination methods except those recommended in this manual, contact Thermo Fisher Scientific to ensure that the proposed method will not damage the device.
28. If a single level of hazard protection is required, make sure to use containers that are tested to Risk Group II or above and that more than one level of protection is provided in the case of materials of a higher Risk Group.
29. Should the unit fail to work to specification, immediately stop centrifuging and allow the rotor to come to a complete stop. Clean and troubleshoot the product according to the instructions under "Trouble Shooting" before any further use of the product. Contact an authorized dealer or the manufacturer if trouble persists.
30. Always work in a manner which endangers neither the user nor any other person.
31. Repairs are to be performed by trained and authorized service personnel only. Contact Thermo Fisher Scientific.
32. Opening the instrument housing or improper use of the instrument voids the warranty. If there is a failure during the warranty period, contact Thermo Fisher Scientific for warranty service.
33. Do not use units that have not been correctly installed or repaired.
34. Do not attempt to disassemble or modify this product.
35. Store at room temperature in a dry area. Do not expose to sunlight, moisture, or extreme temperatures for prolonged periods of time.

WARNING

The laboratory procedure should ensure that no person or hazardous substances are present within a 12" (30 cm) zone around the centrifuge when operating.

WARNING

If this product is not used as recommended by the manufacturer, the overall safety will be impaired.

WARNING

Check MSDS, wear required Personal Protective Equipment, and observe all applicable local and national regulations before dispensing and disposing of hazardous materials.

WARNING

This product does not contain bio-seals as per IEC/EN/CSA 61010-2-020 and cannot provide any level of containment in case of a spill or release of toxic, radioactive, or pathogenic micro-organisms. Thus these materials are not recommended to be used with this product.

Specifications

Environmental Ratings

Note: for Indoor use

Altitude: up to 2000 m

Temperature: 2 °C to 40 °C

Maximum Relative Humidity: 80% for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 40 °C

Technical Features

Capacity: Max. liquid density: 1.2g/cm³ or 1.2g/mL)

Standard round 6-position tube rotor: 12 mL Max

Strip tube (2 x 8-position) rotor: 8 mL Max

Dimensions: 153 mm x 128 mm x 104 mm (LxWxH)

Weight: 1.6 lb (0.73 kg) (with rotor)

1.9 lb (0.86 kg) (with power adapter)

Maximum Speed: 6,000 rpm Approx.(with fully loaded rotor)

Maximum RCF: 2,000 x g

Electrical

Power Adapter Rating

Input: 100-240 V AC, 50/60 Hz, 1.0 A

Output: 12 V DC, 2.2 A

mySPIN 6

Input: 12 V DC

Current (under load): Approx. 0.6 A

Power: Approx. 7.3 W

Accessories:

- Six tube adapters for 0.2 mL tubes and six tube adapters for 0.5 mL tubes; PN: 50146145
- One SnapSpin™ strip tube rotor (can be used in place of installed round rotor); PN: 75004063

NOTE

All accessories are rated for maximum speed.

This product has been tested to the requirements of CAN/CSA C22.2 No. 61010-1 second edition, including Amendment 1, or a later version of the same standard incorporating the same level of testing requirements.

This mark is the confirmation that the unit conforms to the EU guidelines and has been tested according to the following EU Declaration of Conformity.

710213

Conforms to:
EN/IEC 61010-1 2010 ed.3
EN/IEC 61010-2-020 ed.3, 2016
ANSI/UL Std. 61010-1 ed.3, 2012
ANSI/UL Std. 61010-2-020 ed.3, 2016
CAN/CSA C22.2 No.61010-1-12 (R2012)
CAN/CSA C22.2 No.61010-2-020 (R2017)

EC Declaration of Conformity

Manufacturer: Heathrow Scientific® LLC, 440 N Fairway Drive, Vernon Hills, IL 60061 USA

European Contact: Emergo, Europe, Princessegracht 20, 2514AP, The Hague, The Netherlands

Product: mySPIN6

Product Model No(s): 120613

This unit has been constructed and conforms to the following:

EMC Standards and Directives:

IEC 61326-1:2013-b (IEC 61326-1:2012)

FCC Part 15, Subpart B: 2015

ICES-003, Issue 5

2014/35/EU (Low Voltage Directive)

2012/19/EU (WEEE)

RoHS Directive (EU) 2012/19

Safety Standards:

EN/IEC 61010-1 2010 ed.3

EN/IEC 61010-2-020 ed.3, 2016


ANSI/UL Std. 61010-1 ed.3, 2012

ANSI/UL Std. 61010-2-020 ed.3, 2016

CAN/CSA C22.2 No.61010-1-12 (R2012)

CAN/CSA C22.2 No.61010-2-020 (R2017)

Signed:

Gary Kamees

Title: Vice President Product Development/Manufacturing

Location: Vernon Hills, IL

Date of Issue: Jan 2025

UKCA Declaration of Conformity

Manufacturer: Heathrow Scientific® LLC, 440 N Fairway Drive, Vernon Hills, IL 60061 USA

UK Responsible Persons: Emergo Consulting (UK) Ltd., Compass House, Viston Park Histon, Cambridge CB24 9BZ

Product: mySPIN6

Product Model No(s): 120613

This unit has been constructed and conforms to the following:

Directives and Regulations:

Electrical Equipment (Safety) Regulations 2016

The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

The Electromagnetic Compatibility Regulations 2016

EMC Standards:

EN 61326-1:2013-b (IEC61326-1:2012)

ICES-003, Issue 6

Safety Standards:

EN/IEC 61010-1 2010 ed.3

EN/IEC 61010-2-020 ed.3, 2016

ANSI/UL Std. 61010-1 ed.3, 2012

ANSI/UL Std. 61010-2-020 ed.3, 2016

CAN/CSA C22.2 No.61010-1-12 (R2012)

CAN/CSA C22.2 No.61010-2-020 (R2017)

Signed:

Gary Kamees

Title: Vice President Product Development/Manufacturing

Location: Vernon Hills, IL

Date of Issue: Jan 2025

Disposal Responsibilities:

The equipment you purchased may contain hazardous substances that could impact the environment. Per regulations on electronic devices in the European Community, you must use the appropriate disposal systems to avoid exposure of these substances to the environment. The disposal systems will reuse or recycle hazardous materials from your equipment responsibly.

The crossed-out wheeled bin symbol invites you to use those systems.

If you need more information on the collection, reuse and recycling systems, please contact your local or regional waste administration. You can also contact us for more information on the environmental performance of our products.

Operation

1. Remove the contents from the package, examining them carefully for breakage, defects, or missing parts. Save the packaging.
2. Place the mySPIN 6 on a sturdy, level surface away from the edge of the counter and other moving equipment. Make sure the power cord is clear of hot surfaces and other hazards.

WARNING


Spin balanced loads only!

Do not place the tubes asymmetrically nor load unequal volume tubes. Improper placement of tubes will lead to insufficient centrifugation and may cause injury or accident.

3. Prior to every use, inspect unit for cracks, missing pieces, and or any damage that would result in an impaired performance and notify Thermo Fisher Scientific BEFORE use.
4. Verify that the power switch is in the OFF position. Plug the power adapter cord into the rear of the unit and the power cord into an approved outlet. Make sure that the power cord connector can be readily grasped at the rear of the device so that it can be removed to easily disconnect it from the centrifuge.

Turn power switch ON. If there is a smooth whirring sound and the unit accelerates with little or no vibration, mySPIN6 is ready to use. If there are loud or unusual sounds, or excessive vibration, DO NOT OPERATE. Contact Thermo Fisher Scientific® Technical Support.

5. On benchtop prepare sample tubes in a rack so that tubes are filled to equal levels. Close tube lids.
6. To open lid, press down on front latch. The lid will open slightly. Grasp lid and tilt lid toward rear of unit.


WARNING


The laboratory procedure should ensure that no person or hazardous substances are present within a 12" (30 cm) zone around the centrifuge when operating.

7. Tubes of equal weight and size should be placed opposite each other. Use sample or water in tube to provide balance. Proper sample balancing illustrated below.
8. Do not force tubes into the rotor. Use adapters and correctly sized rotors for your application.

SnapSpin Circular Rotor

2 Tubes

3 Tubes








4 Tubes

6 Tubes

SnapSpin Strip Tube


9. Close lid securely. A properly closed lid will grasp the front silicon tab keeping the lid closed. If a clicking sound is heard while the unit is operating, the lid is not fully closed. Do not attempt to operate with the lid removed or open. Turn power switch ON. Unit will not operate with lid open. If there is a smooth whirring sound and the unit accelerates with little or no vibration, mySPIN 6 is ready to use. If there are loud or unusual sounds, or excessive vibration, DO NOT OPERATE. Contact the Thermo Fisher Scientific Technical Support.
10. To stop centrifugation, press power switch and wait until the rotor has completely stopped rotating before opening the lid.

WARNING


Do not open lid until rotor has completely stopped! Avoid severe personal injury or property damage from moving parts.

Rotor Replacement

To remove SnapSpin™ rotor: Press center button firmly and at the same time lift rotor upwards.

To install SnapSpin™ rotor: Set rotor onto the hex shaft. Press center button firmly and downward. Rotor should latch with a click sound.

WARNING

Do not use damaged or cracked rotors!

Maintenance

mySPIN 6 is normally maintenance-free. Clean unit only when it is not plugged into an electrical outlet. When necessary, the housing and rotor can be wiped using a damp cloth and a mild, non-corrosive detergent.

WARNING!

Do not use any solvent on the unit that may attack plastic or cause cracks in the rotors or the unit. Avoid exposure to UV radiation, either natural or man-made. In addition to cracks; discoloration, softening, or crazing of the plastic materials can occur. Corrosion can also occur to the metal parts which are part of the SnapSpin™ rotor feature. This damage can weaken them. Ensure the rotor is protected from corrosion and mechanical damage. The rotor must be cleaned with a neutral cleaning liquid.

NOTICE

Ensure all parts are thoroughly dry prior to operation.

Trouble Shooting

Trouble	Possible Cause	Solution
Not working No light	• Cord Loose • Front ON/OFF switch • Lid not closed	• Check Cord • Press front ON/OFF switch • Close Lid
Excessive vibration	• Rotor not installed correctly • Rotor not balanced correctly	• Reinstall rotor • Check for tubes of equal weight and location (see Operation)
Noise clicking sound	• Lid not fully closed	• Close Lid engaging silicon tab with cover
Rotor falls off	• Rotor not attached	• Reattach rotor (See Rotor Replacement)

Return for Repair

Important

Transporting hazardous materials without a permit is a violation of federal law.

Thermo Fisher Scientific will not accept any mySPIN 6 that is not appropriately cleaned and decontaminated.

In the unlikely event of repair, or when damage to the unit necessitates return, contact Thermo Fisher Scientific and obtain return authorization before sending your mySPIN 6 for service.

Warranty

IN NO EVENT WILL THERMO FISHER SCIENTIFIC OBLIGATION UNDER THIS WARRANTY EXCEED THE PRICE OF THE PRODUCT.

Limited Warranty

Thermo Fisher Scientific warrants that mySPIN 6 will be free from defects in workmanship and material for 3 years from the date of purchase.

If you believe that there is a defect in the product, you must, during the warranty period, notify Thermo Fisher Scientific, provide proof of purchase, and return the product to Thermo Fisher Scientific with a Return Authorization form.

If Thermo Fisher Scientific is properly notified and, after inspection, confirms that there is a defect and the warranty period has not expired, Thermo Fisher Scientific will repair, modify, or replace the product, at its sole option, at no charge.

OTHER THAN THIS LIMITED WARRANTY, THERMO FISHER SCIENTIFIC MAKES NO WARRANTIES, EXPRESS OR IMPLIED, WITH RESPECT TO THE QUALITY OR PERFORMANCE OF THE PRODUCT, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR PARTICULAR PURPOSE WHICH ARE HEREBY DISCLAIMED AND EXCLUDED. THERMO FISHER SCIENTIFIC WILL IN NO EVENT BE LIABLE FOR ANY LOSS OF USE, LOSS OF PROFITS, CONSEQUENTIAL, SPECIAL, EXEMPLARY OR PUNITIVE DAMAGES.

THIS WARRANTY DOES NOT COVER:

ACCESSORY PARTS, SUCH AS RUBBER PARTS, THAT ARE DAMAGED BY LIQUIDS OR MISUSE.

ANY DEFECT OR DAMAGE CAUSED BY IMPROPER OR UNREASONABLE USE OF THE PRODUCT. (mySPIN 6 IS DESIGNED FOR USE ONLY BY TRAINED LABORATORY TECHNICIANS. USE BY ANYONE ELSE WILL VOID THIS WARRANTY.)

ANY PRODUCT THAT HAS BEEN, IN THERMO FISHER SCIENTIFIC SOLE JUDGEMENT, TAMPERED WITH, ALTERED, OR REPAIRED BY ANYONE OTHER THAN THERMO FISHER SCIENTIFIC.

ANY PRODUCT THAT IS INOPERATIVE BECAUSE OF: (a) WEAR OCCASIONED BY USE, (b) NEGLIGENCE, (c) ACCIDENT, (d) INCORRECT MAINTENANCE, OR (e) USE UNDER ABNORMAL CONDITIONS OF TEMPERATURE, DIRT OR CORROSION, OR USE WITH ABRASIVE OR CORROSIVE MATERIALS.

Register your warranty online at www.thermoscientific.com/warranty

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thermoscientific.com/centrifuge

© 2024 Thermo Fisher Scientific Inc. All rights reserved.

Heathrow Scientific is a registered trademark of Heathrow Scientific. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Shown pictures within the manual are examples and may differ considering the set parameters and language.

United States/Canada +1 866 984 3766

Latin America +1 866 984 3766

Austria +43 1 801 40 0

Belgium +32 53 73 42 41

France +33 2 2803 2180

Germany national toll free 0800 1 536 376

Germany international +49 61 84 90 6000

Italy +39 02 95059 552

Netherlands +31 76 579 55 55

Nordic/Baltic States/CIS +358 10 329 220

Russia +7 812 703 42 15

Spain/Portugal +34 93 223 09 18

Switzerland +41 44 454 12 12

UK/Ireland +44 870 609 9203

India +91 22 6716 2200

China +800 810 5118 or +400 650 5118

Japan +81 3 5826 1616

Other Asian Countries +852 2885 4613

Australia +61 39757 4300

New Zealand +64 9 980 6700

Countries not listed +49 6184 90 6000

Thermo
SCIENTIFIC
A Thermo Fisher Scientific Brand