
Leica M651

User manual / Bedienungsanleitung
Mode d'emploi / Manual de empleo

10 705 115 – Version A

Leica
MICROSYSTEMS

Contents

Adjustments before the operation page 7

I. Using the floor stand and the wall- and ceiling mounts


1. Securing the swingarm
2. Brake knobs
3. Positions of tiltable joint and microscope carrier
4. Counterbalancing
5. Positioning the stand at the operating table

II. Using the microscope

6. Positioning the surgical operating microscope
7. Interpupillary distance
8. Eye contact
9. Depth of field
10. Dioptric correction
- 10.1 Dioptric correction without graticule
- 10.2 Dioptric correction with graticule
11. Brake of focusing drive
12. Focusing
13. Illumination
- 13.1 Fitting the bulb
- 13.2 Quick-change lamp mount
- 13.3 Adjusting the control unit
- 13.4 Adjusting the transformer
14. Sterile covers

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Adjustments before the operation

1. Position the stand

⚠ Before attempting to conduct the rollable floor stand to the operating table, secure the swingarm as follows:

- Pull out the knob (21) and engage the safety hook (22) on the pin (12).
- Release the articulation brake knobs (90, 91, 92). Fold up the swingarm and retighten the articulation brake knobs.
- Depress the metal lever at the side of each round foot pedal (13) to release the footbrakes. Now the floor stand can be brought to the operating table, using the handles (17).

On arrival at the chosen location, press hard on the two footpedals (16) to immobilize the base.

When positioning the ceiling mount, use the handswitch to set the motor-driven jib to the working height.

2. Arrange the starting position



Leica stands can be equipped individually for any discipline.

- Make sure that the swingarm is blocked (fig. 2). Slacken the articulation brake knobs (90, 91, 92). Extend the swingarm and retighten the articulation brake knobs (27, 28).
- Suspend first the tiltable joint, and then the microscope carrier, in the perforated plate position selected.
- Retighten the brake knobs.
- Using the drive knobs (52, 53), tilt and incline the surgical operating microscope precisely.

3. Attach the accessories

Leica surgical operating microscopes can be used for a wide variety of applications.

- Ensure that the swingarm is blocked (fig. 2).
- Fit the accessories required (for assistant, photography, TV etc).
- Attach the appropriate binocular tube (a)
- Attach the objective (b) and the wide-field eyepieces for spectacle wearers (c), adding spacing rings and format-indicating graticule if necessary.

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





4. Counterbalance the swingarm

The microscope outfit at the end of the swingarm must be counterbalanced so that the swingarm will remain motionless at any set position.

- Free the safety hook (22) and engage it on the pin (21).
- Steadying the surgical operating microscope by holding the focusing drive (40) in one hand, rotate the tension adjustment knob (25) until the swingarm floats freely without drifting upwards or downwards.

5. Individual adjustments

- Use the brake ring (41) to adjust the free running of the focusing drive.
- Use the brake knobs (90, 91, 92, page 6) to adjust the freedom of lateral movement. If the tiltable joint and the microscope carrier are not provided with drive, engage the articulation brake knobs 93 and 94 also.
- Use the knob (55) to set the interpupillary distance, and set the dioptic positions with the eyelens mounts (63) for the surgeon and, where appropriate, for the assistant.
- Fold back the flexible eyecups if required.
- Adjust the depth of field using the lever (56) on the binocular tube or on the beam splitter.


6. Inspect the illumination

- Connect the power cable.
- Turn on the power switch (70).
- Use the switch (71) to turn on the internal illumination.
- If necessary, change the bulbs in the guide rail of the lamp mount (86).
- If a bulb should fail during the operation, change it by means of the quick-change lamp mount (87).
- Insert filters into the filterholders (43) if appropriate.

7. Final preparations

- Release the safety hook so that the swingarm can move freely (fig. 1).
- Bring the focusing drive (40) to its mid-position (junction of arrows).
- Turn the magnification changer (42) to step 6.4.
- Treat the eyepieces with an anti-misting agent.
- Fit the sterile covers.

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I. Manipulations

1. Securing the swingarm

⚠ Secure the swingarm with the safety hook (22) under the following circumstances:

- During assembly, when fitting the swingarm to the column, and when adding the tiltable joint and the microscope carrier.
- When the accessories are changed
- When the positions of the tiltable joint and the microscope carrier in the perforated plates are changed.
- When the rollable floor stand MS-C is moved
- During work breaks.

The safety hook is engaged by pulling out the knob (21) and slightly moving the swingarm until it clicks into position.


To release the safety hook, push its lever upwards (fig. 2) and slightly move the swingarm until the retaining mechanism controlled by knob (21) clicks into position.

2. Brake knobs


The ease of movement of the brake knobs for the articulation and of the focusing drive handles is continuously adjustable over a wide range. The brake knobs for the articulation (90, 91, 92) are round and black, to distinguish them from the round orange knobs (27, 28) used to disconnect hooked components, and from the cylindrical black handles used to move the instrument around.

- (90) Brake knob for lateral movement around column
- (91) Brake knob for lateral movement of whole swingarm
- (92) Brake knob for rotational movement of front end of swingarm


Adjust the freedom of movement of the brake knobs **before** the operation, in order that these non-sterile knobs may be covered.

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
1. Using the floor stand and the wall- and ceiling mounts

3. Positioning the tiltable joint and the microscope carrier

The tiltable joint and the microscope carrier can be initially positioned as desired with the aid of the system of perforated discs (fig. 1). Each perforated disc offers five possible positions, so that an appropriate combination of positions for the tiltable joint and the microscope carrier can always be found.

Warning:
Never change the configuration whilst the equipment is positioned over the field of operation!

To change the configuration, engage the safety hook (22) (see section 1), and tighten the brakes. Slacken the orange knob (27 or 28) of the swingarm or the tiltable joint with one hand until the end of the steel axle is seen, supporting the microscope with the other (fig. 1). Engage the steel pins in the required position on the perforated disc, push them fully in, and tighten the orange knob (27 or 28).

4. Counterbalancing


The wall and floor stands, and the ceiling mount, must be counterbalanced when they are installed and when their loading is altered by the addition or removal of accessories.

Warning:
Never change the counterbalance whilst the equipment is positioned over the field of operation!

During the counterbalancing procedure, the safety hook (22) must be disengaged (see section 1).

Steady the microscope by means of the positioning handles so that it does not drift up or down (fig. 2).

Rotate the tension adjustment knob until this tendency to drift is eliminated and the swingarm is in equilibrium: a swingarm position will be found at which the knob can easily be turned. The scale shows tension values (low loading, low figure).

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5. Alignment of stand to operating table

The choice of a position for the stand of the surgical operating microscope should take into account the fact that the instrument needs to be readily swung out of the way during phases of the surgical operation when it is not required, and readily recalled on demand.


The location of sterile and non-sterile personnel relative to the sterile controls of the surgical operating microscope must also be considered.

The maximum flexibility of the articulation is best exploited if the swingarm does not point upwards at a highly acute angle.


Before moving the **rollable floor stand MS-C**, release the foot brakes by depressing the metal lever at the side of each round foot pedal (13).

⚠ Ensure that the safety hook and the brakes for the articulation are applied (see sections 1 and 2).


The stand should only be pushed around by means of the handles (17) on its column. When the correct position has been established, apply the brakes (16) of the rollable floor stand MS-C firmly so that they are heard to engage. The **MS-F ceiling mount** can be swung through a maximum of 300°. The height of the jib is motor-adjustable using a handswitch (110). The vertical displacement range is 60 cm. The light key of the on/off switch, and the socket for the handswitch (111), are located on the jib.

1

2

3

6. Positioning the surgical operating microscope

The microscope is best positioned over the field of operation by placing both hands on the handles of the optics carrier (fig. 4), and never by holding the binocular tube. The tiltable joint and the microscope carrier enable it to be positioned accurately by tilting and inclining, using the drive handles (52, 53).

If the surgical operating microscope is attached to an adapter with carrier rod, it can be positioned at right-angles to the microscope carrier (fig. 3). This arrangement is advantageous in gynaecology, urology, plastic surgery, and also if the 180° dual stereo attachment is used.

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II. Using the microscope

7. Interpupillary distance

To maintain sterility, the interpupillary distance is set before the surgical operation is commenced. The binocular tubes T have an interpupillary distance adjustment drive with scale (55) from 52 mm to 76 mm.

It is worth drawing up a list of the interpupillary distances of the various people who will be using the surgical operating microscope.

The other types of binocular tube are adjusted by pushing or pulling the tubes closer together or farther apart.

For the **binocular tube with variable angle of observation**, the individual interpupillary distance can be set by hand within the range 52 mm to 76 mm and retained with the fixing screw (54). The value is shown on a scale.

8. Eye contact







The eye of the observer should be at the point where the illuminating ray paths emergent from the eyepiece converge (fig. 4).

In effect, the more it deviates from this position, the more of the peripheral field is lost.

The correct distance from the eye lens (B) to the eye (A) is 22 mm for the eyepieces designed for these surgical operating microscopes (fig. 4). This distance enables spectacle wearers to retain their spectacles during observation.

The flexible eyecups (C) prevent the spectacle lenses from becoming scratched; if the whole working distance is required, they can be rolled back.

If the user wears no spectacles, the degree of eye contact can be changed by inserting spacing rings (D) (fig. 5, 6). The 6 mm and 9 mm spacing rings, which are screwed into position between the eyecup and the eye lens of the eyepiece (fig. 7), can be used in any combination to obtain the most comfortable viewing position.

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
9. Depth of field

The binocular tubes T, the beam splitter and the 180° dual stereo attachment have a built-in double-iris diaphragm (fig. 1), operated with the lever (56), for increasing the depth of field. For sterile work, the lever can be fitted with a sterilizable extension.


10. Dioptic correction

Ametropia is corrected by rotating the calibrated eyecens mounts of the adjustable eyepieces for spectacle wearers. The scale runs from -5 to +5 diopters. A correction may also be necessary for spectacle wearers, particularly if special glasses for operating or for reading are used.

The dioptic correction is affected by accommodation, which itself is influenced by various factors (brightness, time of day, age, fatigue). If a graticule (reticle) is fitted in the eyepiece (see section 20.5), the dioptic correction value must be determined with this in position as it has an effect optically.

- (56) Lever of double-iris diaphragm
- (57) Screw to lock dioptic setting
- (42) Magnification changer
- (40) Focusing drive
- (55) Interpupillary drive

The parfocality is only assured if the dioptic settings have been set precisely. The image will then remain sharp right from the lowest magnification to the highest.

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
10.1 Dioptic correction without graticule

- a) Set the changer to the **highest magnification**. Ensure that the magnification changer is properly engaged. Place a flat test specimen with sharply-defined features (e.g. a cross marked on paper) at the working distance of the objective used and focus on it using the focusing drive.
- b) Set the changer to the **lowest magnification**. Turn each eyelens mount to position +5 without looking into the eyepieces. Focus on the specimen for each eye by turning each eyelens mount in the minus direction until the image is sharp. At this point the eye is relaxed and focused at infinity, which largely ensures fatigue-free working. Secure the diopter setting with the locking screws (67).
- c) Set the changer to the **highest magnification** again.

Inspect the focusing, adjusting it if necessary with the drive handles. The specimen should now remain in focus for all positions of the magnification changer. If it does not, repeat the procedure.

10.2 Dioptic correction with graticule

- a) Fit the graticule into one of the eyepieces (see section 19), and turn the eyelens mount of **this eyepiece** to +5.
- b) Bring the **graticule** into focus by slowly turning the eyelens mount in the minus direction until the two parallel vertical hairlines are visible separately. Secure the diopter setting with the locking screw (67).
- c) Now correct **the other eyepiece** as described in section 10.1.

Finally, check for freedom from parallax: When the head is moved slightly from side to side, there should be no relative movement between the image of the graticule and that of the subject, i.e. the two must be sharply imaged in the same plane.

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

3

11. Brake of focusing drive

The setting of the brake ring of the focusing drive depends on the weight of the accessories fitted to the optics carrier. Set the focusing drive to its mid-position (where the two black arrows touch, see fig. 2).

Turn the brake ring (41) in the appropriate direction until the required ease of movement is attained.

12. Focusing

For the LEICA M651 surgical operating microscope, the subject of the operation is brought into focus by means of the focusing drive handles (fig. 3). The range is 60 mm. The ease of movement is adjusted before the operation in accordance with the weight of the accessories (see section 11).

The diopter corrections must be set before the final focusing procedure is carried out (see section 10).

– Set the focusing drive to its mid position so that the edge of the cover plate is level with the intersection point of the arrows (dotted line, fig. 2).

To obtain an overall view, engage the **lowest magnification** step, ensuring that it engages completely.

– Manoeuvre the **instrument across to the operating field** and lower it towards the subject by means of the swingarm until the subject is recognizable down the microscope. The working distance approximates to the focal length of the objective selected.

– Use the focusing drive to **focus on the subject**. If the diopter correction has been set properly, the subject will remain sharply in focus for all magnification steps.

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13. Illumination

The illumination is regulated by a control unit on the column of the wall or floor stand or of the ceiling mount. For the universal stand and the swinging-arm stand, a 12 V regulating transformer is used.

The surgical operating microscope has a rack (43) which accepts two filter slides in which 32 mm diameter filters (neutral, daylight, green, or UV-barrier) can be inserted. A heat-absorbing filter is built in.

13.1 Fitting the bulb

A lamp housing for centred 12V/50W halogen bulbs is built into the microscope.

(85) Housing cover

(86) Guide rail of bulb mount



Switch off the lamp. Lift the housing cover. Without touching the bulbs, remove the latter from their packing and slide them, pointing upwards, along the guide rails (fig. 2). Do not touch the quartz envelope.

13.2 Quick-change lamp mount

The LEICA M651 surgical operating microscope has a quick-change lamp mount for two bulbs.

If a bulb should fail during the operation, move the bar of the quick-change lamp mount (87) beneath the lamp housing (fig. 3); the second bulb takes the place to the first and the lamp-change indicator (88) sticks out.

Before starting the next surgical operation, replace the defective bulb and push in the lamp-change indicator.

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

1

13.3 Intensity variation with the control unit

The control unit can be turned about the column.

- (70) Power switch with pilot lamp
- (71) Regulating knob
 - Set the voltage selector to the voltage of the power supply (see section 28).
 - Insert the bulb or bulbs (section 13.1).
 - Connect the power cable (orange).
 - Turn on the power switch (70). The green pilot lamp shines.
 - Adjust the brightness continuously with the regulating knob (71).

2

13.4 Intensity variation with the transformer

By turning the regulating knob clockwise, the transformer is switched on and the voltage, which is shown on a meter, is set steplessly.

- (79) Regulating knob
- (80) Voltmeter
 - Set the voltage selector to the voltage of the power supply (see section 32).
 - Insert the bulb or bulbs (section 13.1).
 - Connect the power cable.
 - Turn the regulating knob (79) clockwise. To avoid overloading the bulb, the click stop at 12V is not to be passed (see section 32).

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14. Sterile covers

All of the controls which the surgeon may need to touch during the operation have sterilizable covers.

- (59) Protective glass for objective, gas-sterilizable
- (51) Clip-on handles for all drive knobs, steam-sterilizable
- (52) Transparent knob cover for magnification changer, steam-sterilizable
- (57) Diaphragm lever extension, steam-sterilizable

The clip-on handles engage when pushed over the metal handle axes. To remove them, press the axial button (fig. 2). They are steam-sterilizable.

The transparent cover for the magnification changer knob is merely pressed into position. It is also steam-sterilizable. The magnification can be read through the cover.




The drive knob (55) for adjusting the interpupillary distance on the «1» version binocular tubes is also to be provided with a sterilizable transparent cover.

Because of its high-quality optical coating, the protective glass for the objective is not steam-sterilizable; we recommend that it be gas-sterilized. A metal ring (58) is threaded on to the objective mount and the protective glass is then snapped into position (fig. 3).

A sterilizable extension (57) is also available for the diaphragm lever.

Sterilizable covers permit the binocular tube with variable angle of observation to be tilted under sterile conditions. They fit over the eyepieces (fig. 5).

Disposable covers are also available.

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II. Components

A complete outfit consists of:

- a Stand
- b M651 optics carrier
- c Binocular tube
- d Two wide-field eyepieces for spectacle wearers
- e Objective
- f Positioning handles, steam-sterilizable
- g Two transparent knob covers, steam-sterilizable
- h Control unit or regulating transformer

15. Optics carrier



The magnification changer in the optics carrier of the LEICA M651 surgical operating microscope produces the total magnification steps 6x, 10x, 16x, 25x and 40x.

The low height of the microscope body (only 81 mm) ensures a comfortable body position during the operation, even with objectives of very long focal length or when the beam splitter is used.

The built-in 12 V/50 W illuminator includes a quick-change lamp mount for two halogen bulbs, a built-in heat-absorbing filter, and two interchangeable filter holders.

The optics carrier also has a focusing drive system with two sterilizable positioning handles, and two transparent sterilizable covers for the magnification changer knobs.

The LEICA M651 optics carrier accepts any of the binocular tubes and objectives, and a range of accessories.

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16. Binocular tubes

Low inclined binocular tubes (fig. 1) and **straight binocular tubes** (fig. 2), both with double-iris diaphragm and with drive for adjusting interpupillary distance, are available. A simplified version of each, without diaphragm, and with manual interpupillary distance adjustment, can also be supplied.

The **15° tube wedge** enables the angle of observation to be increased or decreased by 15°.





The **binocular tube with variable angle of observation** (fig. 3) ensures an ergonomically-sound instrument configuration. It provides for fatigue-free body positions with all types of instrument assembly and with any combination of accessories. The angle of observation can be steplessly varied by to 90° in either direction from the vertical. In addition, the binocular tube can be set for either of two viewing heights. The image remains the right way up, irrespective of the tube inclination.

Fitting a binocular tube:

- Place the binocular tube on the optics carrier.
- Low inclined binocular tube: Push the eye-piece tubes apart.
- Secure the binocular tube to the optics carrier using the clamping screw (45).

To set the interpupillary distance, see section 7.

The **rotatable tube adapter** also fits in the same manner between the binocular tube and either the optics carrier or the beam splitter (fig. 5). It enables the binocular tube to be turned laterally by up to 30° when the instrument is in a tilted position.

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17. Objectives

The objectives, of focal lengths 100, 150, 175, 200, 225, 250, 275, 300, 350 and 400 mm respectively, fit all Leica surgical operating microscopes.

Screw the required objective on to the underside of the microscope in an anticlockwise direction (fig. 1). Location of the thread is helped by initially turning the objective clockwise until it is heard to engage.

1

18. Eyepieces

The wide-field eyepieces for spectacle wearers fit into any binocular tube and into the second observer tubes.

- (61) Sleeve for graticule (reticle)
- (62) Graticule (reticle), in mount
- (63) Eye lens
- (64) Spacing ring
- (65) Flexible eyecup

Fitting the eyepieces:

- Remove the dust caps from the eyepiece tubes and fully slacken screws (66) of the eyepieces. Insert the eyepieces into the tubes so that the screws are opposite to the slits of the spring shackles (fig. 3). Secure each eyepiece with the screw (66).
- Screw the spacing ring or rings in between the eyecup and the eye lens (see section 8). If preferred, the eyecups can be rolled back.
- The diopter correction is set on the eye lens mount of the eyepiece (see section 10).


Fitting the graticule for photo- and cine/TV in the eyepiece:

Unscrew the sleeve from the underside of the eyepiece, slide the mount of the graticule (62) over it, align the pattern of the graticule with the white spot on the sleeve, screw the sleeve back into position, and secure the eyepiece in the eyetube.

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III. Accessories

The Leica surgical operating microscopes can be equipped for dual observation, for an assistant, for photomicrography, for cinemicrography and for television. Before adding accessories, engage the safety hook; after adding them, adjust the counterbalance (section 4).

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19. 180° dual stereo attachment


The 180° dual stereo attachment enables two surgeons to work at the surgical operating microscope simultaneously.

The 180° dual stereo attachment is mounted on the LEICA M651 optics carrier like a binocular tube and is secured with the clamping screw (45). The binocular tubes for the two surgeons are then added. The following combinations have been found useful:

- with microscope horizontal:
 - Operator: Low inclined binocular tube T
 - Assistant: Straight binocular tube T
- with microscope steeply inclined:
 - Operator: Low inclined binocular tube T
 - Assistant: Binocular tube with variable angle of observation
- with microscope positioned normally:
 - Both users: Binocular tubes with variable angle of observation.

The **bracket with carrier rod** (fig. 3, see section 6) is recommended in circumstances where the full cantilever of the stand is used and where more space should be made available for the two surgeons.

The dual stereo attachment incorporates a **double-iris diaphragm** and a **beam splitter** with two lateral outlets which can be used for photography, cine and multi-station viewing (see section 20).

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20. Beam splitter

The beam splitter, which has two lateral outlets, is the nucleus of the modular system of second observer tubes, phototubes and cine tubes.

The beam splitter directs 50% of the light to the observer, the remainder being available at the lateral outlets.

The beam splitter is fitted on the optics carrier so that its diaphragm lever (56) points towards the observer, and is secured by a screw (45). If desired, it can be positioned 15 mm higher by inserting an adapter, so that it is not touched when the sterile magnification changer is operated. Place the binocular tube in position. Unscrew the cover plate (95) from the outlet on the side of the beam splitter to be used. Both outlets can be used simultaneously, e.g. for photography while a second observer is looking down the instrument.

The double-iris diaphragm, which is operated with the lever (56), serves to increase the depth of field and to reduce the light intensity. All of these tubes fit to it in the same manner.




20.1 Tube for second observer

The second-observer tube can be rotated about two axes, and consequently positioned at any angle; the image is kept vertical by turning the knurled sleeve (98).

Mounting:

Introduce the dovetail of the tube (96) into the dovetail mount (97) of the beam splitter and secure it.

Insert the eyepiece, secure it, add the spacing ring or rings (see section 8). Adjust the diopter setting for the eye to be used (see section 10). If the second observer is responsible for framing during photography, fit a graticule into the eyepiece (see section 18). If the surgeon and the second observer use the same eyepiece magnifications, they will have the same total magnification and will see the same field of view.

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20.2 Stereo attachment for second observer

The stereo attachment enables an assistant to observe surgical procedures along with the surgeon and, because of the spatial image, to participate.

Carefully position the stereo attachment in the self-centring dovetail mount of the beam splitter and secure it.

Add the straight binocular tube, the 45° inclined binocular tube or the low inclined binocular tube and fit the eyepieces.

The stereo attachment for second observer is rotatable about three axes and can therefore be positioned in relation to the assistant no matter what the instrument configuration. The rotation of the image is corrected with the ring (98).

20.3 Cine / TV tube

The phototube and the cine tube cause the image produced by the objective to be projected directly into the film plane (single-stage imaging), resulting in larger fields of view, better image contrast and excellent depth of field.



The phototube and the cine tube, and their adapters, are parfocally matched to the binocular tubes and the second observer tubes, so that the surgeon and the assistant can see the image even while the photograph is being taken.

Introduce the self-centring mount of the cine/TV tube carefully into the dovetail of the beam splitter and secure it (fig. 2).

Fit the C-mount adapter (101) to the appropriate camera. Mount this combination on the cine/TV tube and secure it. After the screw (105) has been slackened, the camera can be rotated to any of the four basic positions (at 90° to one another) and secured, taking into account the final orientation of the image.

In cameras which have no automatic brightness regulation, the brightness can be altered by means of the diaphragm lever (56) on the beam splitter.

A focusing and framing graticule is to be placed in that eyepiece of the binocular tube which is nearer to the cine/TV tube. It enables the dioptric correction, which determines the sharpness of the image, to be carried out correctly (see section 20.5).

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

20.4 Phototube

The phototube accepts any single-lens reflex camera (fig. 1). The exposure times obtained with the $f = 250$ mm adapter are shorter than with the $f = 350$ mm adapter, but the latter results in higher magnifications. Be aware that the $f = 250$ mm adapter produces a brightness decrease towards the edge of the field.

Mounting:

Introduce the self-centring mount of the phototube carefully into the beam splitter and secure it.

Fit the camera adapter (T-mount) (100) to either the $f = 250$ mm or the $f = 350$ mm adapter and connect the combined unit to the camera.

Fit this assembly to the phototube (fig. 2).

A focusing and framing graticule is to be placed in that eyepiece of the binocular tube which is nearer to the phototube (see section 18). It shows the format limits for the two different adapter focal lengths. The dioptic correction with the graticule must be carried out accurately, so that the image is also sharp in the film plane (see section 10.2).

If daylight colour films are used, a blue filter can be swung in to compensate the colour temperature (see section 13). Colour films for artificial light are however much to be preferred.

20.5 Focusing and framing graticule (reticle) for photo- and cine/TV tubes

For photomicrography and cinematography it is necessary to insert a focusing and framing graticule into one of the eyepieces of the binocular tube (see section 18). It enables the image to be accurately focused in the film plane, and indicates the area of the subject which will appear on the picture (see section 10.2). Accurate focusing and dioptic correction are essential to the sharpness of the image.





The graticule shows the 35 mm format limits for adapter focal lengths 250 mm and 350 mm.

Fig. 3: Phototube $f = 250$ mm:
Format 35 mm

Fig. 4: Phototube $f = 350$ mm:
Format 35 mm

Fig. 5: Cine/TV tube $f = 107$ mm:
Format 7.16 mm x 9.6 mm
(16 mm projector)

Fig. 6: Cine/TV tube $f = 107$ mm:
Format $7\frac{2}{3}$ in, video $\Delta 6.6 \times 8.8$ mm

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
formato $7,16 \times 9,6$ mm (proy. de 16 mm)

20.6 Dual attachment

Connect the dual attachment to the beam splitter (page 38).

TV camera with C-mount

- Equip the TV camera with a TV objective.
- Fit it into the dual attachment.
- Tighten the clamping screw.
- Release the clamping screw and engage the TV camera in one of the 45° click stop positions.
- Tighten the clamping screw.

Photo camera

- Fit the camera adapter to the camera.
- Screw the photo objective to the camera adapter.
- Fit the camera to the dual attachment.
- Tighten the clamping screw.

Adjusting the brightness for TV


The filter changer contains filters with transmittances 30%, 50% and 100% respectively, one of which can be replaced by the 8% filter supplied with the equipment:

- Remove the TV camera.
- Exchange the filter in the TV outlet.

TV- and photo formats with the dual attachment

fig. 2

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
21. Leica photomicrographic systems

The HV or HU phototube, or the trinocular tube, is the nucleus of the Leica modular photomicrographic system. The eyepiece in the phototube magnifies the intermediate image further (two-stage magnification), and higher image scales are thereby attained.

The trinocular phototubes are matched so that the focusing and framing can be carried out down the binocular tube. The binocular tube with variable angle of observation cannot be used with the HV and HU phototubes.

A format-indicating graticule, valid for all formats, is to be introduced into one of the eyepieces of the binocular.

The built-in double-iris diaphragm enables the depth of field to be individually adjusted.

Leica photoautomats meet every flexibility requirement with regard to the saving of time and materials. Magazines and camera backs for all formats and film materials can be easily and rapidly interchanged.

Detailed information is to be found in the instruction booklets for the Leica modular photomicrographic systems.

A phototube can be used not only for still photography, but also to enable a **cine or TV camera** to be easily fitted. There is a choice of two types of linkage: the TV adapter with C-mount, or the objective $f = 50$ mm with C-mount.

The observer can view the subject down the microscope during the process.

*LEICA M651 surgical operating microscope
with MPS 48/52 Photoautomat*

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IV. Assembly

22. Wall mount MS-A, fixed

⚠ The wall mount is carried by a mounting bracket which must be fitted by an experienced tradesman at the correct height on the wall.

- (1) Mounting bracket
- (2) Column
- (3) Control unit
- (20) Swingarm

Assembly:

- Column

Slide the column through the holes (b) and secure it with the screws (5) (fig. 2). The other screws (6) must be accessible so that the swingarm may be mounted above the bracket.

- Control unit

The control unit may be fitted either above the bracket or below it, depending on the circumstances.

Assembly above the bracket (fig. 3):

Screw the four screws (6) into the column so that they no longer project.

Slide the safety ring along the column and secure it with the screw (4).

Slide the control unit along the column (fig. 3).

- Swingarm

Loosen the four screws (6) until they no longer protrude into the bore.

⚠ Tighten the braking knobs and engage the safety hook (section 1.2).





Lower the swingarm into the column (fig. 4). Tighten the screws (6).

Loosen the screws (5) and the safety ring, and lower the column and the swingarm to an acceptable position.

Tighten the screws (5) and secure the safety ring with the screw (4).

If the control unit is to be fitted **beneath the bracket**, the mounting sequence is as follows:

- Place the column in the bracket.
- Lower the swingarm into the column.
- Lower the column and the swingarm far enough for the control unit to be slid on to the column from beneath.
- Secure the safety ring beneath the control unit.
- Connect the control unit (see section 28).

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23. Floor stand MS-C, rollable

The floor stand MS-C consists of a rollable base, a column and the swingarm.

- (15) Base
- (19) Column
- (3) Control unit
- (18) Metal handle axes
- (17) Handles for pushing the floor stand

Assembly:

- Insert the **column**:
Apply the foot brakes (16) by pressing both of the pedals down hard. Remove the dust of the screws (14). Insert the column (fig. 1). A **20 cm column extension** can be introduced between the column and the rollable base if required for special operating techniques.
- Fit each of the two metal handle axes (18) to the column by means of the two hollow screws, and then add the two screw-on handles (17). Rotate the column so that the handles point out sideways relative to the base, secure the column using the screws (14), and fit the two dust plugs.
- Attach the **control unit**:
Slide the control unit, control knobs upwards, down the column (fig. 2).
- **Fit the swingarm.**




⚠ Disconnect the power cable

- Loosen the four screws (8), and open the control unit.
- Support the cover on account of the electrical earthing wire.

⚠ Tighten the braking knobs on the swingarm (90, 91, 92, section 2).

Ensure that the safety hook (22) is engaged (see section 1).

- Lower the swingarm into the column (fig. 4).
- Tighten the four screws (6).
- Close the control unit.
- Connect the control unit (see section 28).

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
1

24. Ceiling mount MS-F

The fitting of the mounting bracket or ceiling plate, the ceiling mount, the swingarm, the carrier with control unit and the electrical connections to the ceiling mount will be organized by the Leica agency.

The fitting of the microscope carrier and of the tiltable joint is described in sections 25-27, and the connection of the control from the swingarm to the control unit in section 28.

Fuses:

230 V

T 6.3 A, 5 x 20 mm (167656)

T 0.05 A, 5 x 20 mm (167660)

T 0.5 A, 5 x 20 mm (167650)

115 V

T 6.3 A, 6.3 x 32 mm (167673)


T 0.05 A, 6.3 x 32 mm (167676)

T 1.0 A, 6.3 x 32 mm (167669)

25. Tiltable joint

A tiltable joint, can be fitted to the swingarm. It increases the overhang and mobility of the instrument, and has a safety mechanism involving a perforated plate and a crane hook.

- Slacken the orange knob (27) until the end of the steel axle is seen.
 - Hook the tiltable joint on to the axle.
 - Engage the pegs of the tiltable joint at an appropriate position in the perforated disc of the swingarm and tighten the orange knob (27) (the alteration of the peg position on the perforated disc is described in section 3).
- The orange knob (28), which also forms part of the tiltable joint, serves to secure the microscope carrier.

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

26. Optics carrier


26.1 Fitting to the microscope carrier

The LEICA M651 optics carrier, if used with the MS-A mount, the MS-C stand or the MS-F mount, is to be attached to the microscope carrier.


Mounting

- Loosen the two screws (46), remove the cover plate (fig. 2), and lower the focusing drive (40) completely.
- Slide the microscope on the peg (47) of the microscope carrier. Hold the microscope so that the four hollow screws (48) can be inserted, and tighten them.
- Replace and secure the cover plate with the two screws (46).
- To attach the microscope carrier, see section 27.

3


4

26.2 Fitting to the bracket with carrier rod

- Position the bracket with carrier rod in relation to the microscope carrier by means of the peg (47) (fig. 4).
- Fit and tighten the four hollow screws (48).
- Secure the adapter (35) to the optics carrier by means of four hollow screws (fig. 5). The 25 mm diameter hole must be vertically positioned.
- Remove the screws from the carrier rod end-stop (49).
- Slide the optics carrier with adapter on to the carrier rod from beneath and secure it with the clamping screw (54).
- Secure the carrier rod end-stop (49) with the screws; it prevents the instrument from sliding off.

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
1

27. Microscope carrier

- Attach the microscope carrier (with microscope, see section 26) either to the tiltable joint (fig. 1), or directly to the swingarm (fig. 2).
- Slacken the orange knob of the swingarm or the knob (27) of the tiltable joint until the end of the steel axle is seen.
- Hook the tiltable joint on the axle, engage the pegs at an appropriate position in the perforated disc, and retighten the orange knob. In combination with one another, the microscope carrier and the tiltable joint ensure the correct position for the user (see section 3).

28. Control unit

The illuminator is to be connected to the control unit if the wall or floor stands or the ceiling mount are used.

Warning:
Before opening the control unit, always disconnect the power cable.

(73) Voltage selector and fuse holder

(74) Socket for power cable

(75) Socket for lamp cable


- Connect the cable from the lamp housing to the socket on the swingarm.
- Connect the cable from the swingarm to the socket «LAMP» (75) on the control unit (fig. 4).
- Press the cable into the clips, ensuring enough latitude for the movement of the instrument (fig. 1).
- Connect the power cable (orange) to the socket (74).
- Insert the bulb (see section 13.1).
- Connect the power cable to the electricity supply.

Changing fuses


If after the control unit has been properly connected an switched on the illuminator does not shine, first inspect the fuses.

- Pull out the power cable
- Open the cover (73)
- Pull out the holders marked with an arrow. You can tell a blown fuse by its melted glass tube.
- Replace the holders complete with new fuses.
- Open the cover (73)
- Turn the drum so that the appropriate figure 115 V (for 100-120 V) or 230 V (for 200-240 V) is visible in the window when the cover is replaced.

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
V. Care of the instrument

1. Never dismantle optical, mechanical or electronic components unless directions for doing so are given in the instruction booklet.
2. Protect the instrument from dust, damp, and acid or alkali vapours.
3. Do not store any chemicals in the neighbourhood of the instrument.
4. If the instrument no longer functions satisfactorily, call in a specialist or your local Leica agent.
5. **Optics:**
Dirty eyepieces and objectives have a very bad influence on the image quality. Dust is best removed with a pneumatic rubber bulb.
Eyepieces and objectives are to be cleaned with a clean soft cloth;
fingerprints and tenacious marks can be removed with a cloth dipped in alcohol or ether (remove the eyecups first).
Objectives are best cleaned from the centre outwards.
Eyepieces, tubes without eyepieces, and optics carriers without tubes should be protected against dust by means of the covers provided.
6. **Stands:**
Stand, when not in use, should be protected by a dust cover.
Do not grease or oil any parts.
Do not substitute electrical plugs.
7. Only authorized service staff may effect repairs and maintenance.

Guarantee

We guarantee the quality of each of our instruments. Our guarantee covers all faults in materials and manufacture. It does not, however, cover damage resulting from careless or improper handling.

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VI. Optical Data

LEICA M651 surgical operating microscope

Eyepiece Okulaire Okular Ocular	Magnification changer position Vergrößerungswechslerstufe Posición del cambiador de aumentos	Objective / Objectif / Objektiv / Objetivo											
		f=100 mm (1.0x)	f=150 mm (0.67x)	f=175 mm (0.6x)	f=200 mm (0.5x)	f=225 mm (0.44x)	f=250 mm (0.4x)	f=275 mm (0.36x)					
10x/21 B	6	Working distance Distance de travail Arbeitsabstand Distancia de trabajo 91 mm	Working distance Distance de travail Arbeitsabstand Distancia de trabajo 140 mm	Working distance Distance de travail Arbeitsabstand Distancia de trabajo 170 mm	Working distance Distance de travail Arbeitsabstand Distancia de trabajo 190 mm	Working distance Distance de travail Arbeitsabstand Distancia de trabajo 215 mm	Working distance Distance de travail Arbeitsabstand Distancia de trabajo 240 mm	Working distance Distance de travail Arbeitsabstand Distancia de trabajo 265 mm	Total magnification Grossissement total Aumento total	Field diameter Champ visuel Gesichtsfelddurchmesser	Total magnification Grossissement total Aumento total	Field diameter Champ visuel Gesichtsfelddurchmesser	Total magnification Grossissement total Aumento total
	10	6.0x	4.0x	3.4x	3.0x	2.7x	2.4x	2.2x	96.3*	mm	96.3*	mm	96.3*
	16	10.0x	6.7x	5.7x	5.0x	4.4x	4.0x	3.6x	57.8	mm	57.8	mm	57.8
	25	16.0x	10.7x	9.1x	8.0x	7.1x	6.4x	5.8x	36.1	mm	36.1	mm	36.1
	40	25.0x	16.7x	14.3x	12.5x	11.1x	10.0x	9.1x	23.1	mm	23.1	mm	23.1
16x/14 B	6	40.0x	26.7x	22.9x	20.0x	17.8x	16.0x	14.5x	14.4	mm	14.4	mm	14.4
	10	9.6x	6.4x	5.5x	4.8x	4.3x	3.8x	3.5x	64.2	mm	64.2	mm	64.2
	16	16.0x	10.7x	9.1x	8.0x	7.1x	6.4x	5.8x	38.5	mm	38.5	mm	38.5
	25	25.6x	17.1x	14.6x	12.8x	11.4x	10.2x	9.3x	24.1	mm	24.1	mm	24.1
	40	40.0x	26.7x	22.9x	20.0x	17.8x	16.0x	14.5x	15.4	mm	15.4	mm	15.4
25x/9.5 B	6	64.0x	42.7x	36.6x	32.0x	28.4x	25.6x	23.3x	9.6	mm	9.6	mm	9.6
	10	15.0x	10.0x	8.6x	7.5x	6.7x	6.0x	5.5x	43.5	mm	43.5	mm	43.5
	16	25.0x	16.7x	14.3x	12.5x	11.1x	10.0x	9.1x	26.1	mm	26.1	mm	26.1
	25	40.0x	26.7x	22.9x	20.0x	17.8x	16.0x	14.5x	16.3	mm	16.3	mm	16.3
	40	62.5x	41.7x	35.7x	31.3x	27.8x	25.0x	22.7x	10.5	mm	10.5	mm	10.5
	40	100.0x	66.7x	57.1x	50.0x	44.4x	40.0x	36.4x	6.5	mm	6.5	mm	6.5

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VI. Optische Daten

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Microscope d'opération
LEICA M651

Operationsmikroskop
LEICA M651

Microscopio de operación
LEICA M651

Eyepiece Oculaire Okular Ocular	Magnification changer position Vergrößerungswechslerstufe Changeur de grossissements en position	Objective / Objectif / Objektiv / Objetivo			
		f=300 mm (0.33x)	f=350 mm (0.28x)	f=400 mm (0.25x)	
10x/21 B	6	2.0x	1.7x	1.5x	Field diameter Champ visuel Gesichtsfeld Durchmesser mm
	10	3.3x	2.9x	2.5x	Total magnification Grossissement total Aumento total
	16	5.3x	4.6x	4.0x	Field diameter Champ visuel Gesichtsfeld Durchmesser mm
	25	8.3x	7.1x	6.3x	Total magnification Grossissement total Aumento total
	40	13.3x	11.4x	10.0x	Field diameter Champ visuel Gesichtsfeld Durchmesser mm
	6	3.2x	2.7x	2.4x	Total magnification Grossissement total Aumento total
	10	5.3x	4.6x	4.0x	Field diameter Champ visuel Gesichtsfeld Durchmesser mm
	16	8.5x	7.3x	6.4x	Total magnification Grossissement total Aumento total
	25	13.3x	11.4x	10.0x	Field diameter Champ visuel Gesichtsfeld Durchmesser mm
	40	21.3x	18.3x	16.0x	Total magnification Grossissement total Aumento total
16x/14 B	6	5.0x	4.3x	3.8x	Field diameter Champ visuel Gesichtsfeld Durchmesser mm
	10	8.3x	7.1x	6.3x	Total magnification Grossissement total Aumento total
	16	13.3x	11.4x	10.0x	Field diameter Champ visuel Gesichtsfeld Durchmesser mm
	25	20.8x	17.9x	15.6x	Total magnification Grossissement total Aumento total
25x/9.5 B	40	33.3x	28.6x	25.0x	Field diameter Champ visuel Gesichtsfeld Durchmesser mm

* Maximum illuminated field of view:
24 mm for f = 100 mm objective
36 mm for f = 150 mm objective
43 mm for f = 175 mm objective
48 mm for f = 200 mm objective
55 mm for f = 225 mm objective
60 mm for f = 250 mm objective
67 mm for f = 275 mm objective
72 mm for f = 300 mm objective
87 mm for f = 350 mm objective
96 mm for f = 400 mm objective

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VII. Technical data

Only the following accessories may be connected to the control unit:

1. Tested Leica accessories.
2. Accessories passed by an approved testing station as safe to use.

Power supply

Primary
115 V/230 V ~ + 10 %
- 15 %
50 / 60 Hz

Fuses

2 x 1.6 AT (5 x 20 mm)

Safety class

Type B $\wedge$ Class I

Secondary Output

13.5 V = 5 A

Microscope : LEICA M651 for lamp type
12 V / 50 W
Fuse 6.3 AT (5 x 20 mm)

Care!

The control unit was set to **230 V grid voltage** in the factory. To change the grid voltage setting, refer to page 56.

Pull out the cable before you open the control unit.

Use the instrument only in medically-designated rooms in accordance with VDE 0107. Do not use the instrument where there is the danger of explosion.

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Leica Microsystems – the brand for outstanding products

The Leica Microsystems Mission is to be the world's first-choice provider of innovative solutions to our customers' needs for vision, measurement, lithography and analysis of microstructures.

Leica, the leading brand for microscopes and scientific instruments, has grown from five brand names with a long tradition: Wild, Leitz, Reichert, Jung and Cambridge Instruments. Leica symbolizes both tradition and innovation.

Leica Microsystems – an international company with a strong network of customer services.

Australia:	North Ryde BC NSW	Tel. +61 28 17 14 77	Fax +61 28 17 83 58
Austria:	Vienna	Tel. +43 149 54 41 60	Fax +43 149 54 41 630
Canada:	Willowdale/Ontario	Tel. +1 416 497 24 60	Fax +1 416 497 85 16
Denmark:	Herlev	Tel. +45 44 54 01 01	Fax +45 44 54 01 11
Finland:	Espoo	Tel. +35 896 153 555	Fax +35 895 022 398
France:	Rueil-Malmaison Cédex	Tel. +33 147 32 85 85	Fax +33 147 32 85 86
Germany:	Bensheim	Tel. +49 6251 1360	Fax +49 6251 13 61 55
Hong Kong:		Tel. +85 22 56 46 699	Fax +85 22 56 41 63
Italy:	Milan	Tel. +39 02 57 48 61	Fax +39 02 57 40 34 75
Japan:	Tokyo	Tel. +81 35 43 59 600	Fax +81 35 43 59 619
Korea:	Seoul	Tel. +82 23 41 64 400	Fax +82 23 41 64 400
Netherlands:	EK Rijswijk	Tel. +31 70 41 32 100	Fax +31 70 41 32 109
Norway:	Oslo	Tel. +47 67 92 27 00	Fax +47 67 92 27 03
Portugal:	Lisbon	Tel. +35 11 38 89 112	Fax +35 11 38 64 215
Singapore:		Tel. +65 77 97 823	Fax +65 77 30 628
Spain:	Barcelona	Tel. +34 93 49 49 530	Fax +34 93 49 49 532
Sweden:	Sollentuna	Tel. +46 86 25 45 45	Fax +46 86 25 45 10
Switzerland:	Glattbrugg	Tel. +41 1 809 34 34	Fax +41 1 809 34 44
United Kingdom:	Milton Keynes	Tel. +44 1908 246 246	Fax +44 1908 609 992
USA:	Allendale/New Jersey	Tel. +1 20 12 365 900	Fax +1 20 12 365 908

and representatives of Leica in more than 100 countries.

The Business Unit SOM, within Leica Microsystems Ltd, holds the management system certificates for the international standards ISO 9001 and ISO 14001 / EN 46001 relating to quality management, quality assurance and environmental management.

EU-Bevollmächtigter:

Leica Microsystems Holdings Ltd
D-35578 Wetzlar
Tel. +49 6441-290. Fax +49 6441-29 33 99

Leica Microsystems Ltd
Business Unit SOM
CH-9435 Heerbrugg (Switzerland)

Telephone +41 71 727 31 31
Fax +41 71 727 43 43
www.leica-microsystems.com

Leica

M I C R O S Y S T E M S

Microscopes

Compound
Stereo
Surgical
Laser Scanning
Photomicrography
Video Microscopy
Measuring Microscopes

Advanced Systems

Image Analysis
Spectral Photometry
Automated Inspection
Stations
Measurement Systems
Electron Beam Lithography

Laboratory Equipment

Tissue Processors
Embedding Systems
Routine & Immunostaining
Coverslippers
Refractometers

Microtomes

Rotary & Sliding
Cryostats
Ultramicrotomes
EM Sample Preparation