
EN-US Instructions for use PL Instrukcja używania CS Návod k použití HU Használati utasítás SK Návod na použitie RO Instrucțiuni de utilizare HR Upute za upotrebu BG Инструкции за употреба SL Navodila za uporabo SR Uputstvo za upotrebu LV Lietošanas instrukcija ET Kasutusjuhend LT Naudojimo instrukcija RU Инструкция по пользованию MK Упатство за употреба TH คู่มือการใช้งาน MS Arahan penggunaan VI Hướng dẫn sử dụng

WM 68111a 08/2024 EN-US, PL, CS, HU, SK, RO, HR, BG, SL, SR, LV, ET, LT, RU, MK, TH, MS, VI

CE 0197

 **Manufacturer**
Löwenstein Medical
Technology GmbH + Co.KG
Kronsaalweg 40
22525 Hamburg, Germany
T: +49 40 54702-0
F: +49 40 54702-461
www.loewensteinmedical.com

JOYCEclinic Full Face, JOYCEclinic Full Face NV, JOYCEclinic Full Face NV + AAV





Full Face Mask

WM 68111a

LÖWENSTEIN
medical

LÖWENSTEIN
medical

with the warranty periods from date of purchase listed below. The warranty conditions are available on the manufacturer's website. We will also send you the warranty conditions on request.

Please bear in mind that any claim under warranty and liability shall be void if neither the accessories recommended in the instructions for use nor original spare parts are used.

In the event of a claim under warranty, contact your specialist dealer.

Product	Warranty periods
Masks including accessories	6 months

12 Declaration of Conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsaalsweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant regulations of Directive 93/42/EEC governing medical devices. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website.

In the EU: As a user and/or patient, you must report any serious incidents occurring in conjunction with the product to the manufacturer and to the responsible authority.

13 JOYCEclinic ordering information

Article number by number of pieces per packaging unit

Article	1 pc.	10 pcs.
JOYCEclinic Full Face NV Size S Size M Size L	WM 26805 WM 26815 WM 26825	WM 26857 WM 26858 WM 26859
JOYCEclinic Full Face Size S Size M Size L	WM 26840 WM 26850 WM 26860	WM 26866 WM 26867 WM 26868
JOYCEclinic Full Face NV + AAV Size S Size M Size L	WM 26845 WM 26855 WM 26865	WM 26869 WM 26870 WM 26871

5 pcs. per packaging unit

Article	Article no.
JOYCEclinic headgear	WM 26835

1 Operation

Use the fitting template to select the correct mask size (see Figure 4).

The figures show the following steps for operating the mask:

1 Putting on the mask

2 Adjusting the mask

3 Removing the mask

For blind and partially-sighted users

An electronic version of the instructions for use is also available on the manufacturer's website.

2 Introduction

2.1 Intended use

The JOYCEclinic Full Face / JOYCEclinic Full Face NV / JOYCEclinic Full Face NV + AAV full face masks are used as a link between the patient and the ventilator for non-invasive positive pressure respiratory therapy in hospitals/clinics. They are intended for a single patient with respiratory insufficiency (> 30 kg) and spontaneous breathing who may be subjected to non-invasive ventilation with pressure support. The masks may only be used in combination with ventilators that have appropriate alarms and safety systems for possible device failure. Only use JOYCEclinic Full Face NV in combination with ventilators that have an active exhalation system.

2.2 Contraindications

The mask must not be used for patients weighing < 30 kg.

The mask may not be used in the following situations: Necessity for immediate intubation, loss of consciousness, acute vomiting.

The mask may be used in the following situations only with particular caution: pressure points and acute injuries to the facial skin, skin allergies in the facial area, facial or nasopharyngeal deformities, acute pain in the facial area, cough reflex restricted or absent, claustrophobia, acute nausea.

Mask is not suitable for use with a nebulizer or as an anesthetic mask.

If you are not certain whether one of these situations applies to you, consult your healthcare professional. Please observe the contraindications in the instructions for use of your device.

2.3 Side effects

The following side effects may occur with use of the mask: Nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing.

If these side effects occur, contact your healthcare professional.

2.4 Clinical benefit

Transfers the therapeutic efficacy of the therapy device to the patient.

3 Safety

⚠ WARNING! Unusually hazardous situation. If you ignore the following instructions, severe, irreversible or fatal injuries may result.

Risk of injury due to damaged mask parts or those under strain!

- ⇒ Perform a visual inspection before every use and after every cleaning operation. Replace mask parts if necessary.
- ⇒ Note useful life (see section entitled "Technical data").

Risk of injury from re-inhaling CO₂!

- ⇒ Only use the mask when therapy is in progress.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Patients unable to remove the mask themselves must be monitored by a nurse.
- ⇒ Check before every use that the openings of the anti-asphyxia valve are clear.
- ⇒ Do not close off exhalation systems.

Risk of injury due to the use of oxygen!

Oxygen can accumulate in clothing, bed linen, and hair. Supplying oxygen without a safety device can lead to fire.

- ⇒ Use an oxygen safety valve.
- ⇒ Follow the instructions for use for the oxygen supply system.
- ⇒ Do not smoke.
- ⇒ Avoid naked flames.
- ⇒ Ventilate the room well.
- ⇒ Keep mask free from oil and grease.
- ⇒ Observe the safety instructions in the instructions for use of the device.

Risk of injury due to insufficient supply to the patient!

- ⇒ Activate low pressure/leakage alarms on device.
- ⇒ Use the correct mask size and check that it fits tightly.
- ⇒ Monitor patients with restricted spontaneous breathing.

Risk of injury from inadequate cleaning!

- ⇒ Clean mask parts before using for the first time (see section entitled "Cleaning and reprocessing").
- ⇒ Clean the mask regularly.

⚠ CAUTION! Hazardous situation. If you ignore these instructions, mild or moderate injuries may result.

Risk of injury due to allergies!

- ⇒ Be aware of the materials used in the mask (see section entitled "Materials"). Only use a mask following consultation with a healthcare professional.
- ⇒ When selecting the detergent, consider potential allergies.

4 Product description

A diagram of the individual parts can be found on the title page.

1	Headgear	7	Mask body
2	Forehead cushion	8	Mask cushion
3	Forehead support	9	Headgear clip
4	Retaining ring	10	Quick-release cord
5	Elbow	11	Ball joint adapter (for JOYCEclinic Full Face only, not shown)
6	Anti-asphyxia valve		

Compatible devices

In some device combinations, the actual pressure does not correspond to the therapy pressure displayed by the device. Have the device adjusted by a healthcare professional so that the actual pressure in the mask corresponds to the therapy pressure. This adjustment should be made with the type of mask that will be used during therapy.

Exhalation system

Masks with an integrated exhalation system have a gap through which exhaled air escapes.

Only use masks without an integrated exhalation system ("NV", blue retaining ring and elbow) with devices that have an active exhalation system as well as alarms and safety systems for the event of the device failing. If using external exhalation systems, follow the associated instructions for use.

Anti-asphyxia valve (AAV)

If the device fails, the anti-asphyxia valve opens and the patient breathes ambient air.

Quick-release cord

The quick-release cord allows the mask to be released quickly and easily in emergency situations (see SOS illustration).

5 Cleaning and reprocessing

The mask is intended for a single patient and a maximum useful life of 7 days. A change of patient is not permitted.

Clean the mask with an isopropyl alcohol wipe (70 % by vol.) before using it for the first time and if it is lightly soiled.

Note: The mask cannot be safely disinfected using an isopropyl alcohol wipe (70 % v/v). The step described only cleans the mask, it does not disinfect it.

5.1 Disposal

Dispose of the mask in domestic waste. In the clinical environment: Dispose of the mask in accordance with hospital regulations.

6 Troubleshooting

Fault	Cause	Action
Pain due to pressure on the face	Mask too tight.	Adjust headgear so that it is looser.
Draft in the eye	Full face mask fits too loosely.	Adjust headgear so that it is tighter.
	Mask does not fit.	Contact your specialist dealer.
Therapy pressure is not reached.	Mask is not correctly adjusted.	Re-adjust mask.
	Mask cushion is damaged.	Replace mask cushion.
	Circuit is damaged.	Check circuit and the correct fit of circuit.
Therapy pressure not reached.	Anti-asphyxia valve defective.	Replace anti-asphyxia valve.

7 Technical data

	Vented	NV
Dimensions in mm (H x W x D)		
Size S	100 x 145 x 110	100 x 145 x 110
Size M	100 x 155 x 110	100 x 155 x 110
Size L	100 x 165 x 110	100 x 165 x 110
Weight		
Size S	110 g	110 g
Size M	113 g	113 g
Size L	117 g	117 g
Dead space volume		
Size S	225 ml	225 ml
Size M	255 ml	255 ml
Size L	305 ml	305 ml
Tube connection: Tapered connection to EN ISO 5356-1	Ø 22 mm (female)	Ø 22 mm (female)
Flow resistance at 50 l/min at 100 l/min	0.04 hPa 0.19 hPa	0.04 hPa (NV + AAV: 0.05 hPa) 0.11 hPa (NV + AAV: 0.19 hPa)
AAV flow resistance		
Insp. at 50 l/min	0.8 hPa	-
Exsp. at 50 l/min	0.5 hPa	-
Tolerance: ±0.2 hPa		

	Vented	NV
AAV switching pressure		
Open	0.8 hPa	-
Close	3.2 hPa	-

Product class in accordance with directive 93/42/EEC	Ila
Service life	5 years
Useful life	Up to 7 days ¹
Therapy pressure	4 hPa - 50 hPa
Quoted two-figure noise emission value according to ISO 4871:	
Sound pressure level	23 dB(A)
Sound power level	31 dB(A)
Uncertainty factor	3 dB(A)
Temperature:	
Operation	+5 °C to +40 °C
Transport and storage	-20 °C to +70 °C
Standards applied	EN ISO 17510: 2020


¹ Useful life depends on cleaning and on the cleaning agents used, on the amount of time worn daily, on therapy pressure, and on individual secretion of sweat. The use of the second mask cushion (optionally included) will not extend the useful life.

8 Materials

No parts of the mask contain latex, PVC (polyvinyl chloride) or DEHP (diethylhexyl phthalate).

Headgear	CO (cotton), EL (elastane), PA (polyamide), P (polyester), PU (polyurethane)
Headgear clip	POM (polyoxymethylene)
Quick-release cord	PET (polyethylene terephthalate fiber), PA (polyamide)
Quick-release cord clip	POM (polyoxymethylene)
Forehead cushion	SI (silicone)
Retaining ring	PP (polypropylene)
Mask body	HTC (high-temperature copolymer)
Mask cushion	SI (silicone)
Elbow	PC (polycarbonate)
Anti-asphyxia valve	TPE (thermoplastic elastomer), PP (polypropylene)
Rotating joint adapter	PC (polycarbonate)

9 Characteristic pressure/flow curve

10 Markings and symbols

The following markings and symbols may be applied to the device, accessories or packaging.

Symbol	Description
	Manufacturer and, if necessary, date of manufacture
	Unique device identifier (unique product code for medical devices)
	Order number
	Indicates the product is a medical device
	Follow the instructions for use
	Permitted temperature range for transport and storage
	Use by date
	Keep out of sunlight
	Lot number
	Do not reuse
	CE symbol (confirms that the product conforms to the applicable European directives/regulations)

11 Warranty

Löwenstein Medical Technology gives the purchaser of a new original Löwenstein Medical Technology product and of a spare part fitted by Löwenstein Medical Technology a limited manufacturer warranty in accordance with the warranty conditions applicable to the product in question and in accordance

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