



UDR 800

DIGITALIZED REFRACTOR



Brilliantly Conceived & Beautifully Crafted

UDR 800

DIGITAL REFRACTOR FEATURES

Eco-friendly & Design using antibacterial material

Reduced Body Size	Faster Restarting Time	Faster Speed of Lens Change	Reduced Operating Noise
----------------------------------	---------------------------------------	--	--

8.0 inch Full-colored Touch LCD Panel

- It could serve clear distinct UI so that it may serve high resolution graphic.
- Realize UI with various images and distinguish the condition of Operation Panel and data by serving various colors. (Locking system function, Prismoft, Autooft, Autooft function, External data representation and Touch masking function, etc.)

Minimized number of buttons



- 26-Button and simply touch to operate easily.

Short cut key button for often used functions

- Operating test data through jog Dial. (SPH, CYL, AXIS, ADD, VA, Prism, etc.)
- Menu selection through jog Dial.
- Additional functions and data value selection through jog Dial.
- Shortcut key button function.

Perfect Interface with Ref/Keratometer, Chart Projector, Refractor, LCD projector and Lens Meter

- It could serve receiving Data from URK-700 & URK-800 and saving DB.
- It could be connected with A, B, C and D type of ACP-700 including ULC-800 (UDR-700 could be connected only with A type of ACP-700)

Monocular PD measurement
Worldwide Various Chart Support
LCD Monitor Tilting function and Screen reversal function
Embedded printer
Help and Real-time Guide function
It could serve upgrading program and A/S functions by USB port

Simple GUI for User Convenience (Various Vision Test Function Offer)



unicos

UNICOR, Co. Ltd.
 Unicos B/D, #114-1 Munhwa-dong, Jung-gu, Daejeon, Korea
 Tel: +82-42-925-8045 Fax: +82-42-925-8041
 www.e-unicos.com E-mail: marketing@e-unicos.com



Networking System

Efficient data transmission among a series of UNICOR Brands.

Measuring Range

Spherical Lens	-29.00 ~ +28.75D (For general test) +19.00 ~ +18.75D (For uncor cylinder & Prism test) (Step: 0.125D / 0.25D / 0.375D / 0.50D / 0.625D / 0.75D)
Cylinder Lens	0.00 ~ -8.75D (Step: 0.25D / 0.50D / 1.0D / 2.0D / 3.0D)
Cylinder Axis	0 ~ 180° (Step: 1.5 / 5 / 30 / 45°)
Pupil Distance	Far 48 ~ 80 mm (Step: 0.5 / 1.0 mm) Near 45 ~ 75 mm (Step: 0.5 / 1.0 mm)
Working Distance	05 ~ 70 cm (Step: 5 cm)
Rotary Prism	0 ~ 20 Δ (Step: 0.1 Δ / 0.2 Δ / 0.5 Δ / 1 Δ / 2 Δ)
Cross Cylinder	Jackson Cross Cylinder ± 0.25D Jackson Cross Cylinder ± 0.50D Dual Cross Cylinder
Retinoscopy Lens	+1.5D, +2.0D (Test Distance: 67cm, 50cm)

Specifications

Body	362 (W) × 82 (D) × 299 (H), 3.5Kg
Controller	215 (W) × 230 (D) × 226 (H), 1.5Kg
Junction Box	286 (W) × 60 (D) × 239 (H), 1.0Kg
Power Supply	AC 100 ~ 240V, 50/60Hz
Power Consumption	90VA

Auxiliary Lens

Open / Close Lens	
Pin Hole Lens	Φ 1mm
Modiox Rod	Right Eye (Horizontal red), Left Eye (Vertical red)
Red / Green Filter	Right Eye (Red), Left Eye (Green)
Polarized Light Filter	Right Eye (135°, 45°), Left Eye (45°, 135°)
Separating Prism	Right Eye (H&V), Left Eye (H&V) Left Eye (10Δ, BI soon to be added to 0~6Δ)

PD Test Lens	
Fixed Cross Cylinder Lens	Jackson Cross Cylinder ± 0.05D, Axis fixation 90°
Vision Degrees	32°

Designs and details can be changed without prior notice for its improvements.

Distributed by