





The rotary evaporator RV 3 is the ideal entry-level model of the IKA rotary evaporator portfolio. It finds a multitude of uses in the chemical, pharmaceutical, and biotechnology industries, in research and development, in manufacturing and quality assurance, in laboratories, and in plant construction.

- heating bath with 4 I volume for up to 99 °C
- water heating bath with digital temperature display and carrying handles
- mechanical lift end-point safety stop
- locking mechanism: red indicator shows unlocked position of the vapor tube

www.ika.com





- manual lift for precise positioning of the glassware
- adjustable immersion angle
- single-handed manual lift handling, suitable for left and right-handed operators
- stepless speed setting with dial control and speed display
- speed range: 20-300rpm
- low device voltage (24V) ensures user safety
- flask clamping mechanism with integrated push-off function for easy exchange of evaporation flasks
- compatible with the entire range of IKA RV 10 glassware

Scope of delivery •RV 3 ECO •HB eco S099 •RV 05.10 Balljoint clamp •RV 10.88 Clamp NS 29 •RV 10.70 Vapor tube (NS 29/32) •RV 10.8001 Seal





Technical Data

20 - 300	
scale	
±1	
±1	
manual	
150	
440 x 530 x 330	
12	
5 - 40	
80	
IP 20	
100 - 240	
50/60	
75	
	scale ±1 ±1 150 440 x 530 x 330 12 5 - 40 80 IP 20 100 - 240 50/60

