





The next generation of corded small bone power tools is here

# Time after time

For decades, our customers
have relied on Stryker
power tools to allow
surgeons to work at their
best. With our newest
corded small bone power
tool offering, we're
continuing that history of
trust.

**1996** TPS

2003 CORE System **2009** RemB

**2017** CORE2

**2020** TPX











# Introducing TPX

It's everything you've come to expect from Stryker.


New cord and connector

New device coating

Personalized I.D. Touch settings

Improved hand switch





A **new cord** that weighs 25% less and is smaller in diameter, reduced by 10%

**Distinct connection points** enable a confident connection

A stainless steel run/safe hand switch selfaligns and clicks into place, ensuring reliability

New handpiece coating is tolerable to SPD processing. The device is built to be temporarily submerged during cleaning, prior to sterilization

I.D. Touch customizes how the sagittal saw and micro drill responds, giving you the performance and feel you want

## **Versatility** of use

XL oscillating saw NEW

A new mounted blade method provides a more purposeful connection

Universal collet NEW

Reliably grips wires over the product's life with a new, innovative grip design

Hardened stainless steel jaws minimize wear over time

Chuck-style jaws that aid in wire/ pin insertion by self-centering TPX is compatible with a number of attachments, making it suitable for most small bone procedures.

#### **Universal driver**

Control for multiple specialties, all in one driver.



# Compatible with CORE an

# with CORE and CORE 2

However your rooms are set up, TPX meets you where you are.

A surgeon's touch is one of their greatest assets, and our industry-exclusive I.D. Touch software aims to support this. Within the CORE 2 Console power source, I.D. Touch software enables the sagittal saw and micro drill power tool's performance to be adjusted from 0–100%. This allows surgeons to customize how their electric power tool responds when pressure is applied during use—from low to high sensitivity—redefining performance to their liking.





## Maximize the value

## of your investment

### **Customized financing**

Through Stryker's Flex Financial program, we provide a range of smart alternatives designed to fit your organization's needs. Choose from numerous payment structures that can be customized to meet your budgetary needs and help build long-term financial stability.

### Services & support

Stryker's ProCare solutions are designed to decrease downtime, improve employee satisfaction, and optimize the performance of your power tools system. When we focus on what we do best–technology and equipment–your staff can spend more time caring for patients.

# TPX ordering information

### **TPX** accessories

Part number	Description
7400-004-000	TPX cord
7400-009-000	TPX handswitch
1322-000-000	Oscillating saw adaptor wrench
5400-052-000US	CORE 2 Console
5100-007-000	Uni-directional footswitch
5100-008-000	Bi-directional footswitch
5100-015-251	Irrigation clip (for straight and angled
	micro drill attachment)
5100-015-271	Irrigation clip (for long straight and angled micro drill attachment)
5400-031-500	Micro oscillating saw irrigation clip
5400-034-500	Micro sagittal saw irrigation clip
5400-037-500	Micro reciprocating saw irrigation clip
5400-050-001	Single handpiece irrigation cassette
5400-050-002	Dual handpiece irrigation cassette

### **TPX** handpieces

Part number	Description
7400-015-000	.TPX micro drill
7400-031-000	.TPX oscillating saw
7400-131-000	.TPX XL oscillating sav
7400-034-000	.TPX sagittal saw
7400-037-000	.TPX reciprocating saw
7400-062-000	.TPX wire driver
7400-099-000	.TPX universal driver

### **TPX** sterilization solutions

1900-375-010Hybrid insert tray (3/4 size x 6")
6400-277-010Small bone insert tray (3/4 size x 6")
6400-454-000Small bone system sterilization case bundle (3/4 size x 6"), 7102-454-040
case bottom and lid, 6400-277-010 small bone insert tray
6400-554-000Small bone system perforated container bundle (3/4 size x 6"), 7102-550-
030  lid, $7102-552-020  perforated bottom$ , $6400-277-010  small bone insert$
tray
6400-654-000Small Bone System solid container bundle (3/4 size x 6"), 7102-550-030 lid,
7102-652-020 bottom, 6400-277-010 small bone insert tray
1900-454-000Hybrid sterilization case bundle (3/4 size x 6"), 7102-454-040 Case bottom
and Lid, 1900-375-010 Insert tray
1900-554-000Hybrid perforated container bundle (3/4 size x $6$ "), 7102-550-030 Lid,
7102-552-020 Perforated bottom, 1900-375-010 Hybrid insert tray
1900-654-000Hybrid solid container bundle (3/4 size x 6"), 7102-550-030 Lid, 7102-652-
020 Solid bottom, 1900-375-010 Hybrid insert tray
7102-454-040Case bottom and lid
7102-550-030Perforated container lid
7102-552-020Perforated container bottom
7102-652-020Solid container bottom

### **TPX attachments**

Micro drill attachi	nents (7400-015-000)
---------------------	----------------------

5100-015-250	.MD-series medium straight
5100-015-252	.MD-series medium angled
5100-015-270	.MD-series long straight
5100-015-272	.MD-series long angled
5100-015-043	.4:1 straight reducer

#### Universal driver attachments (7400-099-000)

011101010101010101001000000000000000000
4100-110-000AO small drill
4100-131-0001/4" drill chuck (with key)
4100-132-0005/32" drill chuck
4100-135-000Hudson modified trinkle drill
4100-160-000Trinkle drill
4505-133-0001/4" (6.4mm) locking keyless chuck
4505-134-0005/32" (4.0mm) locking keyless chuck
4100-120-000T-latch
4100-062-000Wire collet (0.7 – 1.8mm diameter)
4100-125-000Pin collet (2.0 – 3.2mm diameter)
4100-126-000Adjustable pin collet (2.0 – 3.2mm diameter)
4100-225-000XL pin collet (3.2 – 4.0mm diameter)
4500-125-000Universal collet
4500-126-000Adjustable universal collet

#### Reamer chucks

4100-210-000	Large AO reamer
4100-231-000	1/4" reamer chuck (with key)
4100-235-000	Hudson modified trinkle reamer
4100-260-000	Trinkle reamer
4100-215-000	DHS/DHC reamer

#### Other attachments

4100-400-000 Sagittal saw attachment (for use with the 2296 and 5400	
micro saw blades)	
4100-600-0001:1 bur attachment	
4100-700-000Cordless driver high speed bur attachment	

#### Chuck keys

<i>1</i> -
4100-131-132Universal key for 4100-131 or 4100-132
1331-001-0091/4" chuck key
1355-005-0005/32" thumb grip key

## *s*tryker

D0000021947 Rev AA Copyright © 2020 Stryker

Stryker Orthopaedic Instruments 1941 Stryker Way Portage, MI 49002 USA Phone: +1 269 385 2600

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: TPX, RemB, CORE. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.