

# V.A.C. FREEDOM™

## T H E R A P Y   S Y S T E M

## A portable system for managing wounds



### Easy to Operate

- LCD with menu interface allows for efficient selection of negative pressure wound therapy
- SENSAT.R.A.C.™ Technology monitors and maintains wound site target pressure
- Smart Alarms help ensure patient safety and assist staff in providing care
- 300mL canister minimizes canister changes and is easily removed and replaced

### Patient Comfort

- Long battery life (~12 hours) and lightweight design enable patients to be mobile and active
- Carrying case allows discreet delivery of therapy on the go
- Adjustable rate of dressing draw down intensity for increased comfort

### Therapy Options

#### Negative Pressure:

- 50-200 mmHg

#### Mode of Operation:

- Continuous and Intermittent

#### Intensity:

- Variable

### Specifications

- V.A.C. Freedom™ Therapy Unit: P/N 320000

#### Dimension:

- 6.5"(16.51cm) W x 3.125"(7.9cm) H x 7.5"(19.05cm) D

#### Weight:

- 3.20 lbs. (1.45 kg)

#### Canister volume:

- 300 mL

#### Battery Life:

- Approximately 12 hours  
Electrical Charger/External Supply Input:  
100-240 VAC, 50/60 Hz, 1 Amp

#### Charger/External Supply Output:

- 18 VDC, 2.2 Amp

#### Maximum Electrical Leakage:

- <100 Microamps

#### IEC Classification:

- Medical Equipment
- Equipment not suitable for use in Presence of flammable anesthetic mixture with air, oxygen, or nitrous oxide

### Storage Conditions

Temperature Range..... -4°F (-20°C) to 140°F (60°C)  
Relative Humidity Range.. 10 to 95% non-condensing

### Operating Conditions

Temperature Range..... 50°F (10°C) to 85°F (30°C)  
Relative Humidity Range.. 30% to 75%

**NOTE:** Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a physician and product instructions for use prior to application. Rx only.